


HM Series – Multipurpose Examination lights

The HM Series of lights offers intense white light in a compact flexible design.

HM FEATURES

- >>> Compact, vented shade minimizes heat output
- >>> Stay cool finger guard for easy adjustment of light
- >>> Swivel joint lamp head allows 359° axial rotation and 90° angular adjustment
- >>> Available in wall, c-clamp, and mobile configurations



TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC
Light intensity: Up to 3,000 Footcandles/32,400 Lux
Light source: Halogen bulb – 12V/20W
Light color temp: 3,000 Kelvin
Color: White – RAL 9002

Configurations:	Wall	C-clamp	Mobile
Weight:	6 lbs	2 lbs	15 lbs

Wall, gooseneck arm	Part No.: HM1000-37W
Wall, articulating arm	Part No.: HM28002W
Mobile light, gooseneck arm	Part No.: HM1006-37W

