

Celestial Star™

The Celestial Star™ is ideal for surgical applications and procedures where maneuverability and compact design are key requirements. These areas include emergency room, trauma rooms, and critical care areas. The light provides excellent illumination, and offers market leading reach and flexibility. In addition the Celestial Star™ minimizes shadows and incorporates a balance arm that provides drift free performance. Available in single & dual ceiling, wall, and mobile configurations.



CELESTIAL STAR™ FEATURES

- >>> Unparalleled flexibility with extensive reach and range of motion
- >>> Virtually maintenance free
- >>> Drift free balance arm design
- >>> 114-inch reach from head to toe
- >>> Patented lamp shade reduces heat output
- >>> (3) 35W long life Halogen bulbs (average over 4,000 hours)
- >>> Removable, sterilizable handle
- >>> Available in single & dual ceiling, wall, and mobile configurations

TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC

Light intensity: Up to 6,000 Footcandles/64,500 Lux

Light source: Halogen bulb – Dichroic 12V/3x35W/10°

Light color temperature: Up to 4,100 Kelvin

Medical standards: FDA Class II

Color: White – RAL 9002

Configurations:	Ceiling	Dual Ceiling	Wall	Mobile
Weight:	52 lbs	104 lbs	38 lbs	54 lbs

Ceiling	Part No.: CS2050C
Dual Ceiling	Part No.: CS2050D
Wall	Part No.: CS2050W
Mobile	Part No.: CS2050M
Extra Sterilizable Handle	Part No.: 18544
Color Corrective filter	Part No.: 18676

