

RLM – Magnifiers for closer inspection

This magnifying light offers durability and performance in an economical package. Ideal for detailed inspection during exams and procedures.

RLM FEATURES

- >>> 5 diopter (2.25X) magnification
- >>> 22W Fluorescent ring light tube
- >>> Optional 5 diopter add-on lens available
- >>> Available in C-clamp, or mobile version

TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC

Light intensity: 1,000 Footcandles/10,800 Lux

Light source: 22W Fluorescent tube

Light color temp: 4,000 Kelvin

Colour: White – RAL 9002

Configurations:	C-clamp	Mobile
Weight:	15 lbs	22 lbs

C-clamp	Part No.: RLM-CL-120V
Mobile	Part No.: RLM-CLT-120V
Add-on lens, option	Part No.: RLM-5XMAG

