


Orion LED™

MEDICAL & INDUSTRIAL LIGHTING



Sturdy lamp delivers uniform, high-quality light output – ideal for examination and task light applications.

 **sunnex**®



Orion LED™


The Orion LED light system can be integrated with OEM equipment for technical products that require high quality light output.



Ideal for examination rooms and hospital rooms, the bright Orion LED focused lamp provides near daylight color rendition with reduced heat.



High quality light output, displaying a uniform light picture as well as a high reflection of an object's true color, color rendering index (CRI) +94 – important in most exam and inspection applications.



Features include:

- LED Technology minimizes heat output
- Intensity: 2,800 to 28,000 Lux
- CRI: 94
- Color Temperature: Normal White – Range: 4000
- Flexible gooseneck arm for maximum flexibility
- Dimming capability standard
- Swivel joint head for flexible light positioning
- LED Bulb Life: 50,000 hours at 35% degradation
- 5-Year Warranty

	Standard	Optional
Power Supply	120/240VAC, 60/50hz	Direct Wiring, Panel Mount
Arm Length	Wall Mount: 27" Gooseneck Arm Wall Mount: 44" Articulated Arm Mobile Base: 37" Arm Floor Base: 37" Arm Ceiling Mount: 44" Arm C-Clamp: 38" Arm	10" Rigid Tube
Lamp Head Diameter	2.5" x 4.85" Long	
Lens/Filter	Clear	
Optics	10°	25°
Output	28,000 Lux (max @ 50 cm)	2,800 Lux
Switch	Push Button on Head	
Dimmer	Continuous Rotational	
Cord	Wall Transformer	
Safety Class (Int'l)	IP65	
Housing Color	White	Black, Silver

MOUNTING OPTIONS



Orion LED™

Explanation to OR Series e.g. OR-120W10L

Series	Mounting	Length of Gooseneck	Color	Beam Angle	Input Voltage
OR-	1 = Panel Mount	00 = None 20 = 20 inch 27 = 27 inch 37 = 37 inch	B = Black W = White S = Silver	10 = 10° 25 = 25°	L = 12-36 DC H = 100-240 AC
	2 = High Sq. Base				
	3 = Mobile Base				
	4 = Wall Mount				
	5 = Ceiling Mount				

Sunnex standard offering is white with 10° Optics, IP65 rating and integrated transformer. Optional configurations are available with longer lead time.

STANDARD						OPTIONAL			
Sunnex No.	Mounting	Arm Length	Voltage	Optics	Color	Sunnex No.	Optics	Color	Panel Mount
OR-120W10L	Panel Mount	20"	12-36V DC	10°	White	OR-120__L	25°	B (=Black), S (=Silver)	-
OR-220W10H	Gooseneck with High Square Base	20"	100-240 AC	10°	White	OR-220__H	25°	B (=Black), S (=Silver)	12-36V DC
OR-227W10H	Gooseneck with High Square Base	27"	100-240 AC	10°	White	OR-227__H	25°	B (=Black), S (=Silver)	12-36V DC
OR-237W10H	Gooseneck with High Square Base	37"	100-240 AC	10°	White	OR-237__H	25°	B (=Black), S (=Silver)	12-36V DC
OR-337W10H	Mobile Stand	37"	100-240 AC	10°	White	OR-337__H	25°	B (=Black), S (=Silver)	N/A
OR-400W10H	Wall Mount (Articulating Arm)	-	100-240 AC	10°	White	OR-400__H	25°	N/A	N/A
OR-500W10H	Ceiling Mount	-	100-240 AC	10°	White	OR-500__H	25°	N/A	N/A
OR-627W10H	Wall Mount (Gooseneck/Rigid Tube)	27"	100-240 AC	10°	White	OR-627__H	25°	B (=Black), S (=Silver)	N/A

