

Celestial Star™ MRI Surgical

The Celestial Star™ MRI brings the performance of a surgical light into the MRI suite. Excellent flexibility and reach combined with high performance and space saving design make the Celestial Star™ MRI ideal for these facilities. Unlimited positioning and market leading vertical and horizontal reach allow for optimal functionality around the bore. “No tool” service minimizes risk, and the compact design limits obstruction of the work area around the MRI equipment.

CELESTIAL STAR™ MRI FEATURES

- >>> Tested for MRI compatibility
- >>> 114-inch reach provides head-to-toe coverage
- >>> Drift free balance arm design
- >>> (3) 35W long life halogen bulbs
- >>> Available in ceiling and mobile configurations

TECHNICAL DATA

Power supply: 120V 60 Hz AC or 230-240V 50 Hz AC
Light intensity: Up to 6,000 Footcandles/64,500 Lux
Light source: Halogen bulb – Dichroic 12V/3x35W/10°
Light color temperature: Up to 4,100 Kelvin
Medical standards: FDA Class II
Color: White – RAL 9002

Configurations:	Ceiling	Mobile
Weight:	52 lbs	54 lbs

Ceiling	Part No.: CS 2050C-MRI
Dual Ceiling	Part No.: CS 2050D-MRI
Mobile	Part No.: CS2050M-MRI

